

"Returning independence and dignity". The Whiz with its lily shape, can be used sitting or standing, indoors or outdoors. Just lift a skirt or unzip a fly and you can wee whenever and wherever you want. Simply hold the Whiz against your body; nature and gravity will do the rest.

CE

CE marked and FDA approved
 Manufactured by Jbol Ltd, Oxford,
 England, UK (www.jbol.co.uk)
 Available on NHS prescription
 Continence Foundation approved
 Further testimonials see:
 www.whizfreedom.com

Distributed by JBOL Limited,
 1 Folly Bridge
 Oxford OX1 4LB, U.K.
 Tel: +44 (0) 1865 240572
 Fax: +44 (0) 1865 246565
 www.jbol.co.uk
 Email: info@jbol.co.uk

Being hydrophobic, the Whiz repels all liquid - with one flick it's clean and dry. Made from medical-grade latex free, anti-bacterial plastic, it is also machine washable, eco-friendly and easy to use. With a Whiz dirty public toilets and squat loos are no longer a problem.



whiz[®] freedom


for when nature calls

Allowing women to wee standing or sitting indoors or outdoors, whenever they choose

www.whizfreedom.com

Worlds first and No.1 hydrophobic, anti-bacterial and eco-friendly urine director

- No need to undress
- Hygienic, medical grade and anti-bacterial
- Discreet, portable and re-usable
- Designed not to flowback, splash or spill
- Used by women all over the world





Patents Granted: GB2362577
GB2396819, GB23855276
UK Reg. Design: 2100623, 2103229
International Patent Application:
01/74275, PCT/GB2003/000472
Further patents pending © JBOL Ltd. (May 2006)

Instructions for Use

Thank you for buying a Whiz® Freedom.

You will be amazed at what difference such a small device can make to your quality of life.



1. Hold the Whiz snugly against the body as shown



2-3. Best positions when practising or using the Whiz in a public toilet



4. Using the Whiz without undressing

Important: The Whiz® must be held correctly. The broad area of the opening must be uppermost and held snugly against the body. Its unique lily shape was tested by over 1,400 women in clinical trials and since used successfully by many thousands of women of all ages, from all walks of life. Held correctly the Whiz will not leak and stays completely dry on the exterior. Remember to keep the nozzle pointed downward and away from the body. Do not share a Whiz with others.

Tips on using the Whiz: Ensure the Whiz is held against the body snugly (it's not designed to enter the body). It won't leak, over fill or splash back. Try your Whiz once in the shower or bath, and you will see how easy it is to use – sitting, standing, hovering or crouching – which ever is most comfortable. To use your Whiz standing at a toilet, face the bowl, hold your aim steady and don't forget to raise the toilet seat!



Hold your Whiz this way up



No need to undress: Using the Whiz with a skirt is simple; just hold the fabric out of the line of fire. Underwear need only be pushed to one side. If you're wearing trousers and the zip is long enough, all you have to do is unzip. Trousers with short zips or buttons may have to be lowered slightly. When outside aim downward, keep your back to the crowd and never wee into the wind!

Attaching a Connector Tube: (an accessory from whizfreedom.com): Useful if you want to wee into a bottle or a catheter bag. Pinch the connector tube slightly and wiggle it one cm onto the nozzle of the Whiz. Give it a shake to be sure it's properly attached. Be sure the bottle is lower than the Whiz and that the tube points downward.

If you are bed bound: Attach a connector tube and either lie on your back and elevate your pelvis with a cushion - or roll onto your side and tuck your upper leg closer to your torso than your bottom leg. Hold the Whiz in place and make sure the Whiz and tube slant downward.

NOTE: These two methods require practice. If possible try using the Whiz while sitting on the edge of the bed or kneeling if appropriate. If you are in a chair or a wheelchair (first apply the brakes), edge forward to the front of the seat with the whiz (and tube if applicable) slanting downward. For further tips and testimonials from users, especially hikers, explorers, disabled and mobility impaired wheelchair users visit whizfreedom.com

Caring for your Whiz: The Whiz is made from medical grade anti-bacterial, liquid repellent material, is non-invasive and is designed to be eco-friendly, allowing for multiple use by a single person. Properly used it remains watertight and dry on the outside. All plastics deteriorate over time, and considering the intimate nature of the Whiz, it's advisable to replace it every year, or after extensive use or if it has lost its shape due to long periods rolled up. **NOTE:** urine is sterile when it exits the body and poses no health hazard to the user.

Cleaning: Your Whiz has a special additive that repels all liquid and makes it essentially self-cleansing. After use give it a flick and it will be clean. Once home, the Whiz should be rinsed in water. Even though the Whiz can be re-used by the same person without washing, we recommend that whenever possible it's rinsed in water. The Whiz can be washed in any washing machine at temperatures not exceeding 40°C.

Storage: The Whiz can be rolled or folded temporarily to carry in a pocket or bag. The flexible material will spring back to its original shape as long as it's not kept distorted for long periods. Always store your Whiz unrolled or it may warp permanently.

Updates and Accessories: Enjoy your Whiz and feel free to give us your feedback. Our web site is updated with latest products and accessories and new developments. www.whizfreedom.com.

Technical: The Whiz is made from a medical-grade, anti-bacterial and water repellent and latex free thermoplastic elastomer. It is CE marked, FDA approved and is available on NHS prescription.

The Whiz® is patent and design right granted and protected. The manufacturer JBOL Ltd, have taken every care and attention to ensure that the Whiz Freedom works as effectively as possible. However, under no circumstances will JBOL Ltd be held liable for any damage, injury or inability to use the Whiz Freedom.